Tabla de Cocciones ceramotion® Me

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Paste Opaque	PO	A-D
Paste Opaque Modifier	POM	gingival, orange
Opaque	O	A-D
Opaque Modifier	OM	gingival, orange
Shoulder	SM	A, B, C, D, white, transpa
Gingival	G	1, 2, 3, 4
Base Dentin	BD	A-D
Base Dentin Modifier	BDM	salmon, caramel, ochre, ivory, lemon, vanilla, brown
Dentin	D	A-D
Dentin Modifier Chroma	DM C	A, B, C, orange
Dentin Modifier Fluo	DM F	cream, yellow, orange
Incisal	I	1, 2, 3
Incisal Opal	IO	1, 2, 3
Incisal Transpa	IT	1, 2, 3
Transpa	T	transpa
Incisal Modifier	IM	opal honey, opal white, opal blue, grey, opal grey
Chroma Concept Paste Opaque	CC PO	1 (bleach), 2, 3, 4
Chroma Concept Opaque	CC O	1 (bleach), 2, 3, 4
Chroma Concept Dentin	CC D	1 (bleach), 2 (bleach), 3, 4, 5, 6, 7, 8, 9, 10, 11, 12
Chroma Concept Incisal	CC I	1 (bleach)
Touch Up Paste Opaque	TU PO	Repair 840 °C
Touch Up Base Dentin	TU BD	light, medium, dark
Touch Up Dentin	TU D	light, medium, dark
Touch Up Incisal	TU I	medium, opal, transpa
Correction	C	transpa
Glaze	GL	transpa
Body Stains	B ST	A, B, C
Stains	ST	1 white, 2 vanilla, 3 yellow, 4 orange, 5 pink, 6 purple, 7 blue, 8 grey, 9 olive green, 10 olive yellow, 11 medium brown, 12 red brown, 13 black
Liquids		Me Standard Modelling Liquid, Modelling Liquid, Modelling Liquid +, Paste Liquid, Powder BOL Liquid, Shoulder Liquid, Stains Liquid, Contrast Marker

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